

## Prepare Samples for Illumina Sequencing after library amplification

Illumina sequencing of the resulting plasmids is highly recommended to ensure the resulting plasmid libraries from this protocol are correct, have a large fraction of sgRNAs that match perfectly to the expected library, and the relative fractions of each sgRNA are tightly distributed. You only need ~100 sequencing reads per sgRNA in the library to be able to validate the library (ex. 20M reads for the entire CRISPRi genome-scale library), so you can include your samples as a <5% spike-in to other sequencing runs.

To sequence your libraries, follow Step 3 of the sequencing sample prep protocol with the following modifications:

- If you are spiking into a sequencing sample with high diversity, such as an RNA-seq experiment, you can omit the Set A/B strategy and just sequence with Set A primer pairs and sequencing primer
- For the PCR, perform one 100uL PCR per library. If you have cloned multiple sublibraries, PCR each individually with a different sample index to allow you to detect any cross-contamination or sample mix-ups (although be careful to not introduce this during the PCR!).
- Substitute the genomic DNA input into the PCR with 100ng library plasmid.
- Only perform 15 cycles of PCR.
- Primers for use with single or paired end sequencing. Compatible with both Hiseq 2500 and Hiseq 4000.

### **Crispri TSS (Index) paired end 5'**

aatgatacggcgaccaccgaGATCTACACGATCGGAAGAGCACACGTCTGAACTCCA  
GTCACGCCAATgcacaaaaggaaactcacct

purple= adapter; blue= illumina primer site for sequencing; black= barcode;  
green= complements sg cassette for enrichment

### **Crispr i TSS Common paired end 3'**

caagcagaagacggcatatcgGATCGACTCGGTGCCACTTTTTTC

purple= illumina adapter; black= complements sg for enrichment

### **a) PCR with 100ng of library and 15 cycles only**

- |                        |           |
|------------------------|-----------|
| • Library              | 100uL Rxn |
| • 5X Phusion HF Buffer | (~100 ng) |
|                        | 20 uL     |

- DMSO (100%) 3 uL
- dNTPs (10 mM) 1 uL
- Primer (Index) 100uM 0.25 uL
- Primer (Common) 100uM 0.25 uL
- HF Phusion 1 uL
- Water Fill to 100uL

#### PCR conditions

- 98°C, 30s
- 15 cycles
  - 98°C, 15s
  - 56°C, 15s
  - 72°C, 15s
- 72°C, 10min
- 7°C hold

*After PCR is complete combine reactions in a 1.5 mL nonstick tube. Proceed with SPRIbead purification (or appropriate column purification)*

#### **b) SPRIbead purification**

*Following the index PCR the library is ~264 bp. For purifying with SPRIbeads you will need 1.5 mL nonstick tubes, DynaMag 2 magnet (life tech), 80% EtOH (make fresh), and Qiagen EB.*

- *Add 0.65X reaction volume SPRI beads (65 uL)*
  - *Pipette mix 15X*
  - *Incubate RT for 5 min*
  - *Place on magnet stand 5 min or until supernatant becomes clear*
    - *>300 bp on beads, < 300 bp in supernatant*
  - *Transfer supernatant to new tube. Keep this fraction as your library is here.*
- *Rebind. Add 1.0X your initial reaction volume SPRI*
  - *Pipette mix 15X*
  - *Incubate at RT for 5 min*
  - *Place on magnet stand for 5 min or until clear*
    - *150-300 bp on beads*
  - *Keep beads. Your library is here.*
  - *Wash beads. Add 1mL 80% EtOH. Incubate at RT 2min. Remove EtOH. Repeat for a total of two washes.*
  - *Air Dry. Beads will change from a glossy texture to matte (non-glossy) and will look dry. This can take from 5 to 15 min.*
  - *Elute 50 uL of Qiagen EB*