

Cloning gRNA into Lentiviral Vectors

Annealing:

23 μ L ddH₂O
 1 μ L 100uM top oligo
 1 μ L 100uM bottom oligo
 25 μ L 2X Annealing Buffer (recipe below)
 50 μ L total reaction volume

Incubate at 95°C for 5min. Let oligos gradually anneal while cooling to RT.
 Make a 1:40 dilution of annealed oligos in ddH₂O.
 Annealed oligos can be stored at -20°C.

Ligation: (can be scaled back for high throughput purposes)

10ng digested vector backbone (digested with BstXI and BpI)
 1 μ L 1:40 diluted annealed oligo
 1 μ L fresh 10X T4 ligase buffer
 1 μ L T4 ligase (could probably use less especially in master mix)
 xx μ L ddH₂O
 10 μ L total reaction volume

Incubate at RT for 1-4hrs.
 Transform into 2.5 μ L DH5 α bacteria

2X Annealing Buffer:

200mM Potassium acetate
 60mM HEPES-KOH pH 7.4
 4mM Magnesium acetate

gRNA seq:	NNNNNNNN.....N (N20)				
Ligates into 5'BstXI-BpI3' digested backbone					
		5' end	gRNA seq	3' end	Full Sequence
	<i>5' forward:</i>	TTG	NNNNNNNN.....N (N20)	GTTTAAGAGC	TTGNNNNNNNN.....N (N20)GTTTAAGAGC
	<i>3' reverse (not reverse complemented):</i>	CTTGTTG	NNNNNNNN.....N (N20)	GTTTAAGAGCTAA	CTTGTTGNNNNNNNN.....N (N20)GTTTAAGAGCTAA
	Order 5' forward and reverse complement of 3' reverse above				